

CB Tissue Lysis Buffer (1X) – 100 ML

Catalog No. L6811

Description

A convenient, ready-to-use tissue lysis buffer that is generally used for western blot, Co-IP, Chromatin IP (ChIP), enzyme activity assays, etc.

Lysis buffer contains 50 mM Hepes, 50 mM NaCl, 1% Triton X-100, 5 mM EDTA, 1.5 mM DTT, 10 mM sodium pyrophosphate, 50 mM sodium fluoride, 1 mM sodium orthovanadate, and freshly add protease inhibitor cocktail.

The buffer is sterilized and ready to use for all human or animal tissues of brain, lung, heart, liver, spleen, kidney, aorta, skeletal muscle, etc.

Storage: at 4°C**Shipping Conditions:** Ice pack**Authorized Uses of *Cell Biologics*' Products**

Our products are not authorized for human use, for in vitro diagnostic procedures or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics* without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.