

**C57BL/6 Mouse Embryonic Tissue Lysate (200 ug)**

Catalog No. pIB631

**Description**

Mouse Embryonic Tissue Lysate is prepared from Day 13-19 embryonic tissue of pathogen-free laboratory mouse using Cell Biologics' Tissue Lysis Buffer. Lysis buffer contains 50 mM HEPES, 50 mM NaCl, 1% Triton X-100, 5 mM EDTA, 1.5 mM DTT, 10 mM sodium pyrophosphate, 50 mM sodium fluoride, 1 mM sodium orthovanadate, and freshly add protease inhibitor cocktail (Catalog No. L6811, Cell Biologics).

The buffer is ready to use for all human or animal tissues of brain, lung, heart, liver, spleen, kidney, aorta, skeletal muscle, lymph nodes etc.

**Storage:**

at -80°C upon arrival. Repeated freezing and thawing of lysate is not recommended.

**Shipping Conditions:**

Ice pack or dry ice.

**Authorized Uses of Cell Biologics' Products**

Our products are not authorized for human use, for in vitro diagnostic procedures or for therapeutic procedures. Transfer or resale of any Cell Biologics' products from the purchaser to other markets, organizations or individuals is prohibited by Cell Biologics without the company's written consent. Cell Biologics' Terms and Conditions must be accepted before submitting an order.